


OEM SPO2 Module Manual

(Version V1.0)



Shanghai Berry Medical Electronic Co. Ltd.

OEM SpO2 Module

Main Features

- ✧ Compact size (40x50x8mm), easy installation, high reliability and accuracy
- ✧ Enhanced signal processing to reduce the effect of random artifacts associated with motion and low perfusion
- ✧ Applicable to adult, pediatric and neonatal patients
- ✧ Powered by single 5V, low power consumption
- ✧ Graph includes plethysmogram, pulse strength bargraph and PR digits
- ✧ Communication protocol compatible with BCI
- ✧ Probe interface compatible with BCI probe
- ✧ Converting between RS232 and TTL level
- ✧ Enable all the values displayed in the evaluation software on personal computer

Functions

- ✧ Parameters: SpO2%, Pulse rate
- ✧ One channel plethysmogram
- ✧ Sensor status detection

1. Function

1.1 Measurement parameters

SpO₂, Pulse rate

1.2 Units

SpO₂: %

Pulse Rate: BPM

1.3 Communication

RS232: 4800/115200 Baud Rate

2. Specification

2.1 Measurement Range

SpO₂: 0 -100%

Pulse Rate: 25 ~ 250 BPM

2.2 Resolution

SpO₂: 1%

Pulse Rate: 1BPM

2.3 Accuracy

SPO₂ Accuracy: $\pm 2\%$ (70-100%)

3% (35-70%)

Pulse Rate $\pm 2\%$ (25-250BPM)

2.4 Physical dimensions (without connector)

Length: 5 cm

Width: 4 cm

Height: 1.5 cm

Weight: 20g

2.5 Connection

J1:

PIN1: TXD

PIN2: RXD

PIN3: Gnd

PIN4: +5V

J3: Probe connection

J3_1: R_LED;

J3_2: IR_LED;

J3_3: GND;

J3_4: Opto-;

J3_5: Opto+;

3. Power

Input Voltage: 3-- 5V

Input Current: <50mA

4. Accessory

Adult finger SpO2 sensor

Adult soft sensor

5. Environment

5.1 Temperature

Working: 0 ~ 45°C

Storage: -20~55°C

5.2 Humidity

Working: 30~95%

Storage: 10~100%

6. Standards

EN865 Calibration specification for pulse oximeters

EN60601-1 Safety for medical equipment

EN60601-1-4 General safety standard

7. Host connection

7.1 Serial communication

Communications: BCI protocol

BCI communication settings

- No parity
- 8 data bits
- 1 stop bit
- 4800/115200 baud rate

7.2 BCI communication protocol

Byte	Bit	Description
1	0~3	Signal strength (0~8)
	4	1=no signal, 0=OK
	5	1=probe unplugged, 0=OK
	6	1=pulse beep
	7	Sync bit=1
2	0~6	Pleth
	7	Sync bit=0
3	0~3	Bargraph
	4	Sensor off=1, 0=OK
	5	Pulse research=1, 0=OK
	6	Bit 6 of Byte 3 is bit 7 of the Pulse Rate
	7	Sync bit
4	0~6	Bits 0-6 of Byte 4 are bits 0-6 of the Pulse Rate
	7	Sync=0
5	0~6	SpO2 0-100%
	7	Sync=0